

BT222 Industrial Lab Oven



JPW's line of industrial lab ovens offer performance at a cost effective price. Air convection keeps the temperature uniform throughout the chamber.

Temperature Range

- Ambient +5°F to 500°F

Applications

- Curing
- Mold preheat
- Drying
- Aging
- Electronics
- Shelf Life

Dimensions

- Internal: 24" w x 24" D x 24" H
- External: 50" w x 35" D x 32" H

Exterior

- Painted Steel Gray

Interior

- Aluminized Steel

Circulation System

- One 10" fans generating 400 cfm
- Direct drive motors, Totally Enclosed Fan Cooled (TEFC)
- Air Convection

Doors

- Single Door
- Positive Latch

Heating System

- 3 KW of heat
- Seamless Tubular elements
- Controlled by Solid state relays

Electrical

- 220/1/60
- High limit for safety

Control Panel

- Main power switch
- Heat on/off switch
- Programmable

Options

- Stand
- Disconnect box
- Shelving (unit comes with one shelf)
- Forced Exhaust
- Recorder
- Stainless Interior

Customize Your Order

Call us with your requirement and we can build the product to your custom specification.

6080 State Route 14 · Trout Run, PA 17771

Phone: **570-995-5025**

Mobile: 570-337-3238 – Fax: 570-995-5026



Shown with optional stainless steel interior